## SIMPLE OPERATIOIN







2. Insert into cuvette



3. Insert into channel and test automatically

## **ADVANTAGES**

## **Flexible and efficient** One-step operation

Single test package

#### **Accurate and reliable**

Multi-methodologies 3-level Calibration System Third-party control from Randox

### **Fully automatic testing**

Intelligent temperature control Closed test without pollution risk



# DETAILS High throughput with rapid tests



**High throughout** 

100 Tests/H

Minimum test time: 2min Emergency priority

Build-in standard curve



#### <u>Technical Parameter</u>

Model HP-AFS/1-PLUS W×D×H 555mm×300mm×365mm

Weight 17kg Power AC110~265V

## TEST MENU

Infection	( Cardiac	⊙Tumor Marker
CRP	H-FABP	PSA
SAA	Муо	fPSA
PCT	cTnl	CEA
Hs-CRP	HCY	AFP
Hs-CRP+CRP	BNP	CA199
<b>Renal Function</b>	NT-ProBNP	CA125 ☆
mAlb	LP(a)	CA153 ☆
Cys-C	Lp-PLA2	(2) Gastric Function
NGAL	CK-MB	PGI
α1-MG	D-Dimer	PGII
β2-MG	HDL-C	Нр
RBP	LDL-C	🥍 Rheumatism
Cre	TG	ASO
© COVID-19	TCHO	RF
lgM / lgG Antibody	♦ Immunoglobulin	Anti-CCP
Antigen	IgA	Liver Function
Microelement	IgE	PA
25-OH-VD	IgM	TRF
Ferritin	IgG	ADA
© Diabetes	Fertility	C3
HbA1c (NGSP)	HCG	C4

☆ Under Developing

17 years in IVD **3** R&D centers

150 patents

100 countries

**62**CE certificates







#### **Technology and Innovation**

- Reagent R&D Center in US;
- Instrument technology from Freiburg Germany;
- 150 international patents and 12 invention patents.

#### **Production capacity**

- I 50000 Tests/day;
- Certified by TUV Rheinland with ISO 13485 and CE marks.

### Facility

- World-class high capacity manufacturing facility (10000 sq.m.) with 1500 sq.m.current GMP plants and strict quality control system.
- A new industrial base was established in China Medical City, Jiangsu (35,000 sq.m.)

## A further step of Immunoassay System

Production base: No.3 Building, Fangyi Tech-park, 3 I 3zhujiang Rd, Shijiazhuang city, P.R. China R&D Center: 3938 Trust Way, Hayward, CA94545, USA

Tel: + 86 311 83832628 Fax: + 86 311 83859087

Service Hotline: 400-0191-606

Email: IMD@hipro.us
Website: www.hipro.us





## **Automated Immunoassay System**

Originate in Silicon Valley

High throughout
Full methodology
Automatic intelligent test





